



In Vitro Metabolism Studies of Synthetic Cannabinoids

New Drugs 2014

University Cattolica del Sacro Cuore

May 15, 2014 Rome

Marilyn A. Huestis, PhD

Chief, Chemistry & Drug Metabolism, IRP

National Institute on Drug Abuse, NIH




12 Synthetic Cannabinoid Classes

- Adamantoylindole
- Benzoylindoles
- Cyclohexylphenols
- Classical dibenzopyrans
- Indazole based
- Naphthoylindoles & WIN 55,212-2
- Naphthoylpyrroles
- Naphthylmethylindoles
- Naphthylmethylindenes
- Phenylacetylindoles
- Quinoliny ester (indole based)
- Tetramethylcyclopropyl ketone indoles
- New synthetic cannabinoids constantly emerging

Analytical Challenges

- Vast array of target analytes
- High potency means low concentrations
- Similar chemical structures makes chromatographic separation difficult
- Common metabolites for different targets
- Changing target analyte availability requires constant method addition & validation
- Lack of information about human urinary metabolites of new cannabimimetics






Metabolite Identification


- To determine best metabolites to identify synthetic cannabinoid intake
- To identify potential active metabolites
- To identify potential toxic metabolites
- Supply single known parent drug for hepatocyte metabolism
- Produce phase 1 & 2 metabolites
- Accurate mass provides elemental composition of molecular ion & fragments

Identifying Human Metabolites with Human Liver Microsomes & Hepatic Cell Cultures & High- Resolution Mass Spectrometry

HLM
Preparation
 $T_{1/2}$
determination



1
Incubation
with human
hepatocytes



2
Analysis by
LC-HRMS






3
Software-
assisted data
analysis



Advantages Of High Resolution MS

- Enhanced specificity
- Simplified sample preparation
- Elemental composition: intact $[MH]^+$ & fragments
- Retrospective analysis is possible





Information-Dependent Acquisition

TOF Survey
MS scan

Product
ion scan
MS/MS

Selection
criteria met?

YES


NO

Select Q1,
Scan TOF

Intensity
> 500 cps


Mass defect filter
(range ± 50 mDa)

Dynamic Background Subtraction



Mass Defect Filtering

- Exact mass of JWH-018 $[MH]^+$ 341.1779 Da
→ mass defect 177.9 mDa
- Mass defect filter: filters substances that fall within a specific window of mass defect



Typical Biotransformations Only Minimally Change Mass Defect


+/- 50 mDa window around the parent drug

	Nominal mass	Mass defect
Oxidation	+ 16 Da	- 5 mDa
Hydration	+ 18 Da	- 16 mDa
Demethylation	- 14 Da	- 16mDa
Reduction	2 Da	+ 16 mDa
Glucuronidation	+ 176 Da	+ 32 mDa
Sulfation	+ 80 Da	- 43 mDa



XLR-11




Indole




Tetramethylcyclopropyl
structure

5-Fluoropentyl chain

XLR-11-2'-carboxylic acid



Metabolic XLR-11 Profile






XLR-11 Metabolism


- Biotransformations: hydroxylation, carboxylation, oxidative defluorination, hemiketal & hemiacetal formation, internal dehydration & glucuronidation
- Defluorination leads to UR-144 metabolites
- As expected, more phase 2 metabolites after 3 h incubation
- Main modifications are terminal position on pentyl chain & tetramethylcyclopropyl ring

● ● ● | AKB-48 Metabolism

- APINACA, N-(1-adamantyl)-1-pentyl-1H-indazole-3-carboxamide)
- Novel synthetic cannabinoid found in herbal smoking blends in Japan in 2012
- Strong binding affinity (2x at CB₁ than CB₂)
- DEA Schedule I drug May 2013
- Goal is to profile metabolites of AKB-48 for detection in forensic & clinical laboratories

AKB-48 Metabolism






5F-PB-22

Indole


Ester

Quinoline

5-Fluoropentyl
chain




1st Metabolism Ester Hydrolysis




Phase II: Glucuronidation

Acyl Glucuronides Reactive & Potentially Toxic




5F-PB-22 Metabolites


Intensity, cps




5F-PB-22 Recommended Targets



M17









M19



M22

STS-135 Metabolites





Determining Primary Human Urinary Metabolite Targets of Synthetic Cannabinoids

- Only interested in synthetic cannabinoids where nothing is known about human metabolism
- To date published targets for XLR-11, AKB-48, RCS-4, RCS-8, PB-22, 5F-PB-22, STS-135
- Ultimate goal is to conduct a controlled administration study of synthetic cannabinoid to document human pharmacology
- Compare others based on CB-1 & CB-2 cannabinoid receptor binding



Synthetic Cannabinoids

- Prosecutors & judges indicate need for human pharmacology data
- Designer drugs are the new face of drug abuse
- Present real challenges for clinicians, treatment providers, emergency departments, police, poison control centers, laboratories, legislators, prosecutors, workplace & anti-doping programs
- Important public health & safety issue



Sandrine Pirard
MD, PhD



Karl Scheidweiler,
PhD



Marta
Concheiro, PhD



Allan
Barnes, BS



Megan
Taylor, BS



Ariane
Wohlfarth, PhD



Sebastien
Anizan, PhD



Natalie
Desrosier, MS



Marisol
Castaneto, MS



Kayla
Ellefson, MS



Thank you for your attention!